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

The screenshot shows the 'Image Info' window in ChemiDoc MP. It contains two sections: 'Acquisition Information' and 'Image Information'. In the 'Acquisition Information' section, the 'Exposure Time (sec)' is 0.084 (Optimal Auto-exposure), which is circled in red. In the 'Image Information' section, the 'Data Range (int)' is 36853 - 64896, also circled in red. Other details include: Imager: ChemiDoc™ MP, Serial Number: KNOWNMP, Software Version: 0.9.0.1448, Application: Colorimetric, Excitation Source: White Epi Illumination, Emission Filter: 590/110 Filter, Flat Field Info: Flat Fielded with [14 FLAT_WHITE...], Bias Info: Bias subtraction done with [D:/Ima...], Filter Info: [CHEZ Lens Distortion].

Acquisition Information	
Imager	ChemiDoc™ MP
Exposure Time (sec)	0.084 (Optimal Auto-exposure)
Serial Number	KNOWNMP
Software Version	0.9.0.1448
Application	Colorimetric
Excitation Source	White Epi Illumination
Emission Filter	590/110 Filter
Flat Field Info	Flat Fielded with [14 FLAT_WHITE...]
Bias Info	Bias subtraction done with [D:/Ima...]
Filter Info	[CHEZ Lens Distortion]

Image Information	
Acquisition Date	7/27/2016 8:31:14 PM
User Name	User1
Image Area (mm)	X: 180.0 Y: 144.1
Pixel Size (µm)	X: 65.3 Y: 65.3
Data Range (int)	36853 - 64896

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The screenshot shows the 'Manual' mode settings in ChemiDoc MP. It includes a 'Set Manual Exposure Time' field set to 0.5 sec and a 'Configure Signal Accumulation Mode' section. The 'Configure Signal Accumulation Mode' section has fields for 'First Image' (60 sec), 'Last Image' (180 sec), and 'Images' (3). Below these fields is a graph showing a linear increase in signal intensity over time, with a checkbox for 'Optimal settings' and a legend indicating 'Higher Resolution Longer Wait Time' and 'Lower Resolution Shorter Wait Time'.

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